A/S VESTFROST

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Checklist
Kontrolseddel

UNPE-925 094-675

Product Number: Produktnummer:	
Serial Number: Seriennummer:	
Charging number Ladingsnummer:	
Production check points Kontrol i produktionen	Description: Beskrivelse:
Body Kabin	No scratches, no dents, no rust/ingen ridser, bufer, rust No corner cracks on plastic frame Ingen revner plastkabinens hjørner Condenser holder on back wall (ingenrevner) Fangegreb af holder til rygdensator No contact between tubes and back cover Ingen kontakt mellem rør og kabinet Putty in the sensor holder (E) i hullet til sensor
Check1 Visual inspection	Smooth opening and closing/Åbning og lukning Hinges (proper fixed) / Fjænger (af hængler) No gap between gasket and frame Ingen åbning mellem tætningslæbe og karm Controller Ingen revner i end-plaster/ingen revner i ender stykker Proper alarm function/ Tilbedingslarm Ingen ribber display til termometer/kontrol
Check2 Electrical inspection	Insulate strength testing is carried out with 1000V 2 sec. Approved with 5mA max. Insulation resistance is carried out at 100% at 500V.c Approved with 20kOhm min.
Check3 Cooling test	Cooling test is carried out 100% and approved with each cooling criteria (Cooling curve available)
Comments: Kommentarer:	
Place of inspection Inspektionssted	Check date: Godkendelses dato: Check by: Godkendt af:
Name: A/S Vestfrost	11/6/21
Address: Falkøvej 12, DK-6705, Esbjerg Ø, Denmark	Sors
Address:	

0191127 Rev 02

VESTFROST SOLUTIONS

VLS 054 SDD

GB Instructions for use
FR Consignes d'utilisation
ES Instrucciones para el uso
RUS Инструкция по эксплуатации
CHI 使用说明

VESTFROST SOLUTIONS

Technical manual VLS 024 / 054 / 094 / 154 SDD

HOW TO MAINTAIN YOUR VACCINE AND WATERPACK FREEZER

- Correct placement of the appliance**
 - Location of the appliance
 - Check for level
 - Check for ventilation
 - Check for power supply
 - Check for door opening
 - Check for door locking
 - Check for door gasket
 - Check for door seal
 - Check for door latch
 - Check for door handle
 - Check for door lock
 - Check for door key
 - Check for door lock cylinder
 - Check for door lock plug
 - Check for door lock pin
 - Check for door lock bolt
 - Check for door lock nut
 - Check for door lock washer
 - Check for door lock spring
 - Check for door lock pinion
 - Check for door lock gear
 - Check for door lock rack
 - Check for door lock pinion gear
 - Check for door lock pinion gear teeth
 - Check for door lock pinion gear backlash
 - Check for door lock pinion gear wear
 - Check for door lock pinion gear damage
 - Check for door lock pinion gear misalignment
 - Check for door lock pinion gear binding
 - Check for door lock pinion gear noise
 - Check for door lock pinion gear vibration
 - Check for door lock pinion gear heat
 - Check for door lock pinion gear oil
 - Check for door lock pinion gear grease
 - Check for door lock pinion gear lubrication
 - Check for door lock pinion gear maintenance
 - Check for door lock pinion gear safety
 - Check for door lock pinion gear security
 - Check for door lock pinion gear reliability
 - Check for door lock pinion gear durability
 - Check for door lock pinion gear performance
 - Check for door lock pinion gear efficiency
 - Check for door lock pinion gear effectiveness
 - Check for door lock pinion gear quality
 - Check for door lock pinion gear quantity
 - Check for door lock pinion gear value
 - Check for door lock pinion gear cost
 - Check for door lock pinion gear benefit
 - Check for door lock pinion gear risk
 - Check for door lock pinion gear return
 - Check for door lock pinion gear investment
 - Check for door lock pinion gear opportunity
 - Check for door lock pinion gear solution
 - Check for door lock pinion gear answer
 - Check for door lock pinion gear response
 - Check for door lock pinion gear action
 - Check for door lock pinion gear reaction
 - Check for door lock pinion gear result
 - Check for door lock pinion gear effect
 - Check for door lock pinion gear impact
 - Check for door lock pinion gear influence
 - Check for door lock pinion gear power
 - Check for door lock pinion gear force
 - Check for door lock pinion gear strength
 - Check for door lock pinion gear energy
 - Check for door lock pinion gear work
 - Check for door lock pinion gear effort
 - Check for door lock pinion gear labor
 - Check for door lock pinion gear time
 - Check for door lock pinion gear money
 - Check for door lock pinion gear value
 - Check for door lock pinion gear cost
 - Check for door lock pinion gear benefit
 - Check for door lock pinion gear risk
 - Check for door lock pinion gear return
 - Check for door lock pinion gear investment
 - Check for door lock pinion gear opportunity
 - Check for door lock pinion gear solution
 - Check for door lock pinion gear answer
 - Check for door lock pinion gear response
 - Check for door lock pinion gear action
 - Check for door lock pinion gear reaction
 - Check for door lock pinion gear result
 - Check for door lock pinion gear effect
 - Check for door lock pinion gear impact
 - Check for door lock pinion gear influence
 - Check for door lock pinion gear power
 - Check for door lock pinion gear force
 - Check for door lock pinion gear strength
 - Check for door lock pinion gear energy
 - Check for door lock pinion gear work
 - Check for door lock pinion gear effort
 - Check for door lock pinion gear labor
 - Check for door lock pinion gear time
 - Check for door lock pinion gear money
- Protect the power supply**
 - Check for power supply
 - Check for power supply voltage
 - Check for power supply current
 - Check for power supply frequency
 - Check for power supply phase
 - Check for power supply ground
 - Check for power supply protection
 - Check for power supply safety
 - Check for power supply security
 - Check for power supply reliability
 - Check for power supply durability
 - Check for power supply performance
 - Check for power supply efficiency
 - Check for power supply effectiveness
 - Check for power supply quality
 - Check for power supply quantity
 - Check for power supply value
 - Check for power supply cost
 - Check for power supply benefit
 - Check for power supply risk
 - Check for power supply return
 - Check for power supply investment
 - Check for power supply opportunity
 - Check for power supply solution
 - Check for power supply answer
 - Check for power supply response
 - Check for power supply action
 - Check for power supply reaction
 - Check for power supply result
 - Check for power supply effect
 - Check for power supply impact
 - Check for power supply influence
 - Check for power supply power
 - Check for power supply force
 - Check for power supply strength
 - Check for power supply energy
 - Check for power supply work
 - Check for power supply effort
 - Check for power supply labor
 - Check for power supply time
 - Check for power supply money
- Loading vaccines**
 - Check for loading vaccines
 - Check for loading vaccines quantity
 - Check for loading vaccines quality
 - Check for loading vaccines value
 - Check for loading vaccines cost
 - Check for loading vaccines benefit
 - Check for loading vaccines risk
 - Check for loading vaccines return
 - Check for loading vaccines investment
 - Check for loading vaccines opportunity
 - Check for loading vaccines solution
 - Check for loading vaccines answer
 - Check for loading vaccines response
 - Check for loading vaccines action
 - Check for loading vaccines reaction
 - Check for loading vaccines result
 - Check for loading vaccines effect
 - Check for loading vaccines impact
 - Check for loading vaccines influence
 - Check for loading vaccines power
 - Check for loading vaccines force
 - Check for loading vaccines strength
 - Check for loading vaccines energy
 - Check for loading vaccines work
 - Check for loading vaccines effort
 - Check for loading vaccines labor
 - Check for loading vaccines time
 - Check for loading vaccines money
- Preventive maintenance**
 - Check for preventive maintenance
 - Check for preventive maintenance schedule
 - Check for preventive maintenance tasks
 - Check for preventive maintenance frequency
 - Check for preventive maintenance duration
 - Check for preventive maintenance location
 - Check for preventive maintenance personnel
 - Check for preventive maintenance equipment
 - Check for preventive maintenance materials
 - Check for preventive maintenance safety
 - Check for preventive maintenance security
 - Check for preventive maintenance reliability
 - Check for preventive maintenance durability
 - Check for preventive maintenance performance
 - Check for preventive maintenance efficiency
 - Check for preventive maintenance effectiveness
 - Check for preventive maintenance quality
 - Check for preventive maintenance quantity
 - Check for preventive maintenance value
 - Check for preventive maintenance cost
 - Check for preventive maintenance benefit
 - Check for preventive maintenance risk
 - Check for preventive maintenance return
 - Check for preventive maintenance investment
 - Check for preventive maintenance opportunity
 - Check for preventive maintenance solution
 - Check for preventive maintenance answer
 - Check for preventive maintenance response
 - Check for preventive maintenance action
 - Check for preventive maintenance reaction
 - Check for preventive maintenance result
 - Check for preventive maintenance effect
 - Check for preventive maintenance impact
 - Check for preventive maintenance influence
 - Check for preventive maintenance power
 - Check for preventive maintenance force
 - Check for preventive maintenance strength
 - Check for preventive maintenance energy
 - Check for preventive maintenance work
 - Check for preventive maintenance effort
 - Check for preventive maintenance labor
 - Check for preventive maintenance time
 - Check for preventive maintenance money
- SCP's - Standard Operating Procedures**
 - Check for SCP's
 - Check for SCP's quantity
 - Check for SCP's quality
 - Check for SCP's value
 - Check for SCP's cost
 - Check for SCP's benefit
 - Check for SCP's risk
 - Check for SCP's return
 - Check for SCP's investment
 - Check for SCP's opportunity
 - Check for SCP's solution
 - Check for SCP's answer
 - Check for SCP's response
 - Check for SCP's action
 - Check for SCP's reaction
 - Check for SCP's result
 - Check for SCP's effect
 - Check for SCP's impact
 - Check for SCP's influence
 - Check for SCP's power
 - Check for SCP's force
 - Check for SCP's strength
 - Check for SCP's energy
 - Check for SCP's work
 - Check for SCP's effort
 - Check for SCP's labor
 - Check for SCP's time
 - Check for SCP's money

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

INSTRUCTION & SERVICE MANUAL

Congratulations:

Now, you are the pride owner of "Bluebird " Automatic Voltage Stabilizer which protects your equipment's from the undesired fluctuation in mains voltage.

To use the stabilizers, some simple points for reference are mentioned below:

- First connect the stabilizer to mains. Please be sure that there is no loose connection and earthing in the mains socket is proper.
- The stabilizer is provided with four Indication lamp

MAINS ON	-BLUE LED
OUTPUT ON	-GREEN LED
HIGH & LOW CUT OFF	-RED LED
- There is delay of 200 second initially whenever, we switch ON the stabilizer & restart the stabilizer
- There is provision of Quick start switch, if we do not want delay then we can push this Quick start switch to ON the stabilizer & its Output.
- Please note that there is One MCB provided at the back of stabilizer. The MCB should be at ON position to start the stabilizer & its Output.
- There is one Schuko/world socket provided at the back of stabilizer to connect your equipment/Refrigerator through it.
- This stabilizer work for the Input voltage variation between 150-280V. The stabilizer will not give Output if Input voltage is lower than 150V-5V (i.e. 145V) & above 280V+ 5V (i.e. 283V-285V). It is likely that voltage remains very low during peak hours at some remote areas. Hence stabilizer may be in Cut off position if the line voltage is below 150V-5V (i.e. 145V). When the Input voltage comes to a safe range, the stabilizer will give Output Automatically.

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

VESTFROST SERVICE APP

1. Download and install app from App store/Play store

PROCOLOS DE INSTALACIÓN

ILR's

SDD's

Annex E – Installation checklist On-grid refrigerator / freezer

Note: The installation technician must fill in this checklist for each completed installation.

Mains powered refrigerator installation checklist		Date:
Country:	Region:	
Installation technician:		
Installation company:		
Address:		
Tel:		
Email:		
Note: All checks must be satisfactory before the installation is handed over to the user.		
CHECK 1 – System description		
1.1	Qualified supplier:	
1.2	Refrigerator Model:	
	Serial number: _____	
	Model reference: _____	
	Product number: _____	
CHECK 2 – Shipment details		
2.1	Was the shipment damaged?	<input type="checkbox"/> Yes <input type="checkbox"/> No
2.2	Were any components missing or under-supplied?	<input type="checkbox"/> Yes <input type="checkbox"/> No
2.3	Have damaged/missing/under-supplied parts been replaced?	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable
	Comments:	
CHECK 3 – Functionality test		
3.1	Have all tests been carried out in accordance with the qualified supplier's instructions?	<input type="checkbox"/> Yes <input type="checkbox"/> No
3.2	30-day temperature logger installed and functional	<input type="checkbox"/> Yes <input type="checkbox"/> No
3.3	Is GSM coverage given at the health facility?	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable
3.4	Are all system components functioning properly?	<input type="checkbox"/> Yes <input type="checkbox"/> No
	Comments:	
CHECK 4 – Training		
4.1	Number of health facility staff trained in usage of refrigerator	
4.2	Number of staff trained in preventive maintenance of refrigerator	

CHECK 5 – Documentation	
5.1	Check if the following documentations are supplied
	- User manual for all system components <input type="checkbox"/> Yes <input type="checkbox"/> No
	- Technician's manual for all system components <input type="checkbox"/> Yes <input type="checkbox"/> No
	- Installation manual <input type="checkbox"/> Yes <input type="checkbox"/> No
	Comments:
CHECK 6 – Overall conclusions and recommendations	
6.1	Recommendation: <input type="checkbox"/> Pass <input type="checkbox"/> Fail
	If FAIL, list outstanding work still required:
	If PASS, the installation can be handed over to the user.
Installation technician signature: _____	
Health center responsible signature: _____	
Name and designation: _____	
Rubber Seal Of the Health facility/ District Health officer	
Date: _____	

Annex E – Installation checklist SDDs

Note: The installation technician must fill in this checklist for each completed installation.

Solar refrigerator installation checklist		Date:
Country:	Region:	
Installation technician:		
Installation company:		
Address:		
Tel:		
Email:		
Note: All checks must be satisfactory before the installation is handed over to the user.		
CHECK 1 – System description		
1.1	Qualified supplier:	
1.2	Refrigerator Model:	
	Serial number: _____	
	Model reference: _____	
	Product number: _____	
1.3	Solar panel:	
	Model Ref: _____	Quantity of panels: _____
CHECK 2 – Shipment details		
2.1	Was the shipment damaged?	<input type="checkbox"/> Yes <input type="checkbox"/> No
2.2	Were any components missing or under-supplied?	<input type="checkbox"/> Yes <input type="checkbox"/> No
2.3	Have damaged/missing/under-supplied parts been replaced?	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable
	Comments:	
CHECK 3 – Solar panel installation		
3.1	Panel Orientation	
	Has the panel been installed at the correct angle towards the equator? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Do shadows fall on the panel between 9:00am and 3:00pm? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	If YES, the system FAILS - if possible remove the shade or the panel must be moved.	
3.2	Panel support structure	
	Are roof fixings in place and are they adequate? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Have theft-deterrent fasteners been used? <input type="checkbox"/> Yes <input type="checkbox"/> No	
3.3	Lightning protection:	
	Has the lightning protection circuit been correctly fitted? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Has the earth electrode been correctly fitted? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Has lightning protection system been tested for electrical continuity? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Comments:	
CHECK 4 – Array cabling and installation		
4.1	Only the solar array cable provided by supplier was used for installation? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Are all electrical connections and array cables concealed and properly protected? <input type="checkbox"/> Yes <input type="checkbox"/> No	
	Comments:	

CHECK 5 – Functionality test	
5.1	Have all tests been carried out in accordance with the qualified supplier's instructions? <input type="checkbox"/> Yes <input type="checkbox"/> No
5.2	30-day temperature logger installed and functional <input type="checkbox"/> Yes <input type="checkbox"/> No
5.3	Is GSM coverage given at the health facility? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not applicable
5.4	Are all system components functioning properly? <input type="checkbox"/> Yes <input type="checkbox"/> No
	Comments:
CHECK 6 – Training	
6.1	Number of health facility staff trained in usage of SDD
6.2	Number of staff trained in preventive maintenance of SDD – periodical cleaning of solar panels etc
CHECK 7 – Documentation	
7.1	Check if the following documentations are supplied
	- User manual for all system components <input type="checkbox"/> Yes <input type="checkbox"/> No
	- Technician's manual for all system components <input type="checkbox"/> Yes <input type="checkbox"/> No
	- Installation manual <input type="checkbox"/> Yes <input type="checkbox"/> No
	Comments:
CHECK 8 – Overall conclusions and recommendations	
8.1	Recommendation: <input type="checkbox"/> Pass <input type="checkbox"/> Fail
	If FAIL, list outstanding work still required:
	If PASS, the installation can be handed over to the user.
Installation technician signature: _____	
Health center responsible signature: _____	
Name and designation: _____	
Rubber Seal Of the Health facility/ District Health officer	
Date: _____	

Protocol de funcionamiento

Functionality test procedure - Vestfrost CCE AC appliances

Check	AC	Outcome
1	Wait for the delay set off voltage stabilizer 3-6 min	Yes No
2	Green power light is ON	Yes No
3	Feel by hand if the compressor is running after 5 minutes of power ON	Yes No
4	Feel by hand if the compressor discharge tube is heating up	Yes No
5	Note temperature on data logger – Place it in basket holder (Picture to be captured in Vestfrost APP)	°C
6	Feel by hand if the inner lining/evaporator plate is cooling after 1 hour	Yes No

Functionality test procedure - Vestfrost CCE SDD appliances


Check	SDD	Outcome
1	Green power light is ON	Yes No
2	Check that controller display is ON	Yes No
3	Feel by hand if the compressor is running after 25 minutes of power ON	Yes No
4	Feel by hand if the compressor discharge tube is heating up	Yes No
5	Note temperature on data logger – Place it in basket holder (Picture to be captured in Vestfrost APP)	°C
6	Feel by hand if the inner lining/evaporator plate is cooling after 1 hour	Yes No

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
Page 2 of 78

GPS Coordinates of installed units

VESTFROST SOLUTIONS

BIOMEDICAL COLD CHAIN MERCHANDISING WINE FOOD & BEV

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
FRIDGE-TAG 2E

FRIDGE-TAG 2E

3 days ago, or on the 19th September 2012, the temperature went below the limit at 18.21h. Due to the “Super-Jump” Function, the Fridge-tag® 2 jumps directly to the latest Alarm event.


The min. temperature was - 1.1°C and the total duration below the limit was 1h 35 min.

After 30 seconds without pressing any button or plug in the USB port the Fridge-tag® 2 gets into LOC mode for 10 minutes to adapt again the storage temperature.

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Read-Out Document of the Q-tag Fridge-tag® 2

Identification number: 1
 Date and time of report creation: 18.12.2012 08:50h
 Upper alarm limit: Above +8.0°C for 1h
 Lower alarm limit: Below +2.0°C for 1h
 Low battery since: 29.11.2012





Nr.	Date (dd MM yyyy)	Events*	Mean temp.	Lower alarm limit		Upper alarm limit		Signature / notes		
				Status	Min. temp.	Duration out of range	Alarm time		Status	Max. temp.
1	Today		+4.9°C	In progress	+4.3°C	0min	In progress	+6.2°C	0min	
2	17.12.2012		+5.3°C	ok	+4.3°C	0min	ok	+7.3°C	0min	
3	18.12.2012		+5.2°C	ok	+4.8°C	0min	ok	+5.7°C	0min	
4	15.12.2012		+5.2°C	ok	+4.8°C	0min	ok	+5.7°C	0min	
5	14.12.2012	t	+7.9°C	ok	+4.2°C	0min	ALARM	+28.4°C	7h 23min	12:12h
6	13.12.2012		+4.8°C	ok	+4.1°C	0min	ok	+6.3°C	0min	
7	12.12.2012	t	+5.0°C	ok	+4.1°C	0min	ok	+9.8°C	20min	
8	11.12.2012		+5.0°C	ok	+4.3°C	0min	ok	+6.3°C	0min	
9	10.12.2012		+4.9°C	ok	+4.3°C	0min	ok	+5.8°C	0min	
10	09.12.2012		+4.9°C	ok	+4.4°C	0min	ok	+5.6°C	0min	
11	08.12.2012		+4.8°C	ok	+4.3°C	0min	ok	+5.5°C	0min	
12	07.12.2012		+4.9°C	ok	+4.3°C	0min	ok	+7.3°C	0min	
13	06.12.2012		+4.8°C	ok	+4.4°C	0min	ok	+6.1°C	0min	
14	05.12.2012		+4.8°C	ok	+4.3°C	0min	ok	+5.7°C	0min	
15	04.12.2012		+5.0°C	ok	+4.4°C	0min	ok	+6.2°C	0min	
16	03.12.2012		+5.0°C	ok	+4.0°C	0min	ok	+18.9°C	50min	
17	02.12.2012		+4.5°C	ok	+4.0°C	0min	ok	+5.3°C	0min	
18	01.12.2012		+4.6°C	ok	+3.9°C	0min	ok	+5.2°C	0min	
19	30.11.2012		+4.7°C	ok	+3.9°C	0min	ok	+5.6°C	0min	
20	29.11.2012		+4.0°C	ALARM	-21.0°C	4h 6min	10:07h	ok	+21.3°C	53min
21	28.11.2012		+5.0°C	ok	+2.9°C	0min	ok	+7.2°C	0min	
22	27.11.2012		+3.2°C	ok	+2.7°C	0min	ok	+3.7°C	0min	
23	26.11.2012		+3.2°C	ok	+2.9°C	0min	ok	+3.7°C	0min	
24	25.11.2012		+3.1°C	ok	+2.9°C	0min	ok	+3.5°C	0min	
25	24.11.2012		+3.1°C	ok	+2.9°C	0min	ok	+3.5°C	0min	
26	23.11.2012	t	+3.4°C	ok	+2.9°C	0min	ok	+20.2°C	40min	
27	22.11.2012		+3.2°C	ok	+2.9°C	0min	ok	+3.8°C	0min	
28	21.11.2012		+3.1°C	ok	+2.9°C	0min	ok	+3.5°C	0min	
29	20.11.2012		+1.8°C	ALARM	-24.2°C	5h 2min	14:35h	ok	+4.0°C	0min
30	19.11.2012		+3.2°C	ok	+2.8°C	0min	ok	+3.8°C	0min	

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