

# MF 214

## VACCINE / ICEPACK FREEZER

The MF series provides vaccine storage and icepack freezing with stable temperature control in the range -15°C to -25°C. The units' thick insulation ensures exceptional hold-over time in case of power cut.

### TrueBlue Series A NEW GENERATION OF PRODUCTS



-15°C  
-25°C

### DIMENSIONS

|                                    |               |
|------------------------------------|---------------|
| Height, mm                         | 910           |
| Weight, mm                         | 1120          |
| Depth, mm with handle              | 700           |
| Gross weight, kg (incl. packaging) | 93            |
| Net weight, kg                     | 71,30         |
| Packed dimensions HxWxD, mm        | 1025x1185x727 |
| Shipping volume m3                 | 0,88          |

### SPECIFICATIONS

|                                         |               |
|-----------------------------------------|---------------|
| Gross volume                            | 171           |
| Icepack freezing capacity, kg/24h       | 7.2           |
| Icepack storage capacity, L             | 180 x 0.6     |
| Temperature range (+43C AMB)            | -15° to -25°C |
| Hold-over time, hours                   | 2.9           |
| Energy consumption - stable, kWh/24h    | 3.00          |
| Energy consumption - cool down, kWh/24h | 2.96          |
| Refrigerant                             | R600a         |
| WHO PQS code                            | E003/025      |
| Climate class                           | T             |

### FEATURES

|                                    |           |
|------------------------------------|-----------|
| Storage baskets, no                | 2         |
| Temperature control                | Automatic |
| Lock + keys                        | YES       |
| Adjustable legs                    | YES, 10cm |
| Junction box                       | YES       |
| Safety thermostat                  | NO        |
| Appliance plugs                    | E, F, G   |
| Industrial grade locking mechanism | YES       |
| Fridge-tag 2E                      | YES       |

### LOADING QUANTITIES

|                              |         |
|------------------------------|---------|
| Qty. per 20' / 40' container | 24 / 60 |
| Qty. per truck (EU standard) | 72      |

### VOLTAGE REGULATION RANGE

|                                                                                        |                  |
|----------------------------------------------------------------------------------------|------------------|
| 220-240V/50/60Hz                                                                       | 110-278V/50/60Hz |
| 100-127V/50/60Hz                                                                       | 82-159V/50/60Hz  |
| Integrated and extended voltage stabilizer as standard in line with PQS E007/VS01-VP.6 |                  |

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

All MF AC models are PQS certified according to the newest performance standards.

WHO PQS (Performance, Quality, and Safety) ensures that immunization equipment meets international standards. This certification guarantees quality, reliability, and safety. Non-PQS certified equipment risks poor performance, safety issues, and higher costs.