



# MF 314

## VACCINE / ICEPACK FREEZER

The MF series provides vaccine storage and icepack freezing with stable temperature control in the range -15°C to -25°C. The units' thick insulation ensures exceptional hold-over time in case of power cut.

Watch product video here: [MF SERIES](#)

All MF AC models are PQS certified according to the newest performance standards.



DIMENSIONS		VALUE
Height, mm		910
Weight, mm		1560
Depth, mm with handle		700
Gross weight, kg (incl. packaging)		117,2
Net weight, kg		80
Packed dimensions HxWxD, mm		1025x1620x727
Shipping volume m3		1,21
SPECIFICATIONS		VALUE
Gross volume		281
Icepack freezing capacity, kg/24h		7,2
Icepack storage capacity, pcs		256 x 0,6L
Temperature range (+43C AMB)		-15° to -25°C
Hold-over time, hours		4,0
Energy consumption - stable, kWh/24h		4,23
Energy consumption - cool down, kWh/24h		4,29
Refrigerant		R600a
WHO PQS code		E003/023
Climate class		T
FEATURES		VALUE
Storage baskets, no		3
Temperature control		Automatic
Lock + keys		YES
Adjustable legs		YES, 10cm
Junction box		YES
Safety thermostat		NO
Appliance plugs		E, F, G
Industrial grade locking mechanism		YES
Fridge-tag 2E		YES
LOADING QUANTITIES		VALUE
Qty. per 20' / 40/ container		18 / 42
Qty. per truck (EU standard)		48
VOLTAGE REGULATION RANGE		VALUE
220-240V/50/60Hz		110-278V/50/60Hz
100-127V/50/60Hz		82-159V/50/60Hz