



VLS 064 RF AC

ICELINED REFRIGERATOR & ICEPACK FREEZER

This model has a built-in heating element ensuring that the appliance can maintain stable temperatures between +2°C to +8°C in all ambient temperatures, setting a new industry standard beyond the PQS requirement.

Watch product video here: [ICELINED COMBI SERIES](#)



WHO PQS (Performance, Quality, and Safety) ensures that immunization equipment meets international standards. This certification guarantees quality, reliability, and safety. Non-PQS certified equipment risks poor performance, safety issues, and higher costs.

| DIMENSIONS | |
|--|------------------|
| Height, mm | 910 |
| Width, mm | 720 |
| Depth, mm with handle | 700 |
| Gross weight, kg (incl. packaging) | 104,5 |
| Net weight, kg | 81,5 |
| Packed dimensions HxWxD, mm | 1025 x 785 x 727 |
| Shipping volume m3 | 0,58 |
| SPECIFICATIONS | |
| Gross volume, L | 75 |
| Net volume, L | 52,5 |
| Temperature range (+43°C AMB) | +2°C to +8°C |
| Icepack freezing capacity, kg/24h | 1,6 |
| Icepack storage capacity, pcs | 6 x 0,6 L |
| Hold-over time, hours (+43°C AMB) | 45 |
| Energy consumption - stable, kWh/24h | 0,63 |
| Refrigerant | R600a |
| Freeze protection, grade | A |
| Climate class | T |
| WHO PQS code | E003/070 |
| FEATURES | |
| Storage baskets, no | 1 |
| Temperature control | Automatic |
| Lock + keys | YES |
| Junction box | YES |
| Adjustable legs | YES, 10cm |
| Safety thermostat | YES |
| Appliance plugs | E, F, G |
| Industrial grade locking mechanism | YES |
| Fridge-tag 2E | YES |
| EMS as option | YES |
| LOADING QUANTITIES | |
| Qty. per 20' / 40' container | 44 / 92 |
| VOLTAGE REGULATION RANGE | |
| 220-240V/50/60Hz | 110-278V/50/60Hz |
| 100-127V/50/60Hz | 82-159V/50/60Hz |
| Integrated and extended voltage stabilizer as standard in line with PQS E007/VS01-VP.5 | |

All our VLS AC models are PQS certified according to the newest performance standards. Our fully integrated remote equipment monitoring system (EMS) can be added as an option for this model.