



VLS 504A AC

ICELINED REFRIGERATOR

The icelined refrigerator series has a built-in heating element ensuring that the appliance can maintain stable temperatures between +2° to +8°C in all ambient temperatures, setting a new industry standard beyond the PQS requirement.

Watch product video here: [ICELINED REFRIGERATOR SERIES](#)

WHO PQS (Performance, Quality, and Safety) ensures that immunization equipment meets international standards. This certification guarantees quality, reliability, and safety. Non-PQS certified equipment risks poor performance, safety issues, and higher costs.



| DIMENSIONS | |
|-----------------------------------------|------------------|
| Height, mm | 910 |
| Width, mm | 1560 |
| Depth, mm with handle | 700 |
| Gross weight, kg (incl. packaging) | 193,7 |
| Net weight, kg | 172 |
| Packed dimensions HxWxD, mm | 1025x1620x727 |
| Shipping volume m3 | 1,21 |
| SPECIFICATIONS | |
| Gross volume, L | 314 |
| Net volume, L | 242 |
| Temperature range (+43°C AMB) | +2°C to +8°C |
| Hold-over time, hours (+43°C AMB) | 55,3 |
| Energy consumption - stable, kWh/24h | 0,64 |
| Energy consumption - cool down, kWh/24h | N/A |
| Refrigerant | R600a |
| Freeze protection, grade | A |
| Climate class | T |
| WHO PQS code | E003/113 |
| FEATURES | |
| Storage baskets, no. top | 6 |
| Storage baskets, no. bottom | 4 |
| Temperature control | Automatic |
| Lock + keys | YES |
| Junction box | YES |
| Adjustable legs | YES, 10cm |
| Safety thermostat | YES |
| Appliance plugs | E, F, G |
| Industrial grade locking mechanism | YES |
| EMS as an option | YES |
| LOADING QUANTITIES | |
| Qty. per 20' / 40/ container | 18 / 42 |
| VOLTAGE REGULATION RANGE | |
| 220-240V/50/60Hz | 110-278V/50/60Hz |
| 100-127V/50/60Hz | 82-159V/50/60Hz |

All our VLS AC models are PQS certified according to the newest performance standards. Our fully integrated remote equipment monitoring system (EMS) can be added as an option for this model.

Updated 10/2024