

## DECLARATION OF CONFORMITY

Document no.: DK250311  
Manufacturer: A/S Vestfrost  
Address: Falkevej 12  
6705 Esbjerg Ø  
Denmark  
Product designation: See Enclose No. 1.

A/S Vestfrost confirm that laboratory fridges and freezers with explosion protection are appropriate for the storage of explosive and highly flammable materials within the meaning of EU Directive:  
**2014/34/EU (ATEX)**

Requirements has been assured by compliance with:

**EN 60079-0 and EN 60079-15**

The system is marked with following:

**Protection for storage: CE  II 3G nA IIB T6**

**A/S VESTFROST**


2025.03.11

Jørgen Dyrholm

---

name

Application Certification Manager  
Esbjerg Denmark



---

signature

Application Certification Manager  
Esbjerg Denmark

Enclosure No. 1.

Declaration of Conformity

Date 11-03-2025

<u>Identification</u>	<u>Model</u>	<u>Type</u>
VESTFROST	R-90i	+2°C to +20°C Refrigerator
VESTFROST	R-138i	+2°C to +20°C Refrigerator
VESTFROST	R-427i	+2°C to +20°C Refrigerator
VESTFROST	F-90i	-5°C to -25°C Freezer
VESTFROST	F-138i	-5°C to -25°C Freezer
VESTFROST	F-427i	-5°C to -25°C Freezer
VESTFROST	VTS098	Premium Ultra Low Temperature -86° Freezer
VESTFROST	R-281 Essential	+2°C to +10°C Refrigerator
VESTFROST	R-350 Essential	+2°C to +10°C Refrigerator
VESTFROST	VT 73 -25°C	-25°C Freezer
VESTFROST	VT 146 -45°C	Low Temperature Freezer -45°C
VESTFROST	VT 206 -45°C	Low Temperature Freezer -45°C
VESTFROST	VT 306 -45°C	Low Temperature Freezer -45°C
VESTFROST	VT 406 -45°C	Low Temperature Freezer -45°C
VESTFROST	VT 546 -45°C	Low Temperature Freezer -45°C
VESTFROST	VT 75 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 147 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 207 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 307 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 307 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 407 -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 547 VF -60°C	Low Temperature Freezer -60°C
VESTFROST	VT 78 -86°C	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	VT 208 -86°C	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	VT 308 -86°C	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	VT 408 -86°C	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF- C74i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-C198i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-C296i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-C383i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-37i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-7i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-51i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-91i	-20°C to -86°C Ultra Low temperature freezer
VESTFROST	ULTF-500	-40°C to -90°C Ultra Low temperature freezer
VESTFROST	ULTF-700	-40°C to -90°C Ultra Low temperature freezer