



**COLD CHAIN
EQUIPMENT**

VLS 096A RF SDD

SOLAR DIRECT DRIVEN TRUE COMBI

The True Combi series can maintain stable temperatures between +2°C to +8°C in all ambient temperatures. The included solar panel kit of 4 x 180W, is connected plug & play, directly to the appliance.

Watch product video here: [TRUE COMBI SERIES](#)



All VLS SDD RF models have been tested and PQS approved according to the newest performance standards. Our fully integrated remote equipment monitoring system (EMS) can be added as an option.

DIMENSIONS		VALUE
Height, mm		940
Weight, mm		1510
Depth, mm with handle		750
Gross weight, kg		190
New weight, kg		150
Shipping volume m3		0.58
REFRIGERATOR SPECIFICATIONS		VALUE
Gross volume, L		147
Net volume, L		110
Temperature range (43C AMB)		+2 to +8C
Autonomy, hours		114
Refrigerant		R600a
Freeze protection, grade		A
Climae class		T
WHO PQS code		E003/132
FREEZER SPECIFICATIONS		VALUE
Gross volume, L		28 x 0.6
Water pack storage capacity, L		2.4
Water pack freezing capacity, L		-18 to -28C
Temperature range (43C AMB)		R290 + R600a
Refrigerant		R600a
FEATURES		VALUE
Storage baskets, no		2
Temperature control		Automatic
Lock + keys		YES
Junction box		YES
Adjustable legs		YES, 10cm
Safety thermostat		YES
EMS as option		YES
Industrial grade locking mechanism		YES
Automatic drainage system		YES